

Single cell RNAseq (10X) sample delivery instructions

1. Cell concentration should be in the range of 700-1200 cell/ul. Please bring at least 30,000cells (recommended 50,000cells).
2. Cell viability should be >90%. For viability of <70%, dead cell removal must be performed prior to bringing the cells.
3. Follow the '[cell prep](#)', '[cell counting](#)' and '[dead cell removal](#)' protocols as needed for more detailed instructions.
4. For customized and optimized dissociation protocols please follow the link <https://www.worthington-biochem.com/tissuedissociation/>